FORTUNA[®] OPTIMAT[®] 2

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F5h Program error

The program in the memory which was to be started has a faulty parameter (volume = 0, for example).

F6h Profile is running

The command can not be run because a profile is still running.

F7h Dialogue is running

The command can not be run, because a dialogue with the user is in progress.

F8h Manual displacement active

The command can not be run since the system is being operated manually.

F9h Wrong flask

The automatic flask sensing system has detected a different type of flask compared to the one transferred. The transferred parameter will thus be ignored.

FAh Unknown system parameter

The transferred system parameter is not used by the system and can thus not be processed.

54h Dispensed volume

Query: 54h
Reply: 54h, xx

xx: Dispensed quantity, format 8 digit BCD, unit µl. The transferred value corresponds to the last value for the volume as indicated by the LCD display. Since for the output the last (not indicated) digit is rounded off, the transmitted value will be too high by 0.5 µl.

55h Conversion factor Steps / ml

Query: 55h
Reply: 55h, xx

xx: Conversion factor, BCD, 8 digits

56h Maximum value ml/s

Query: 56h
Reply: 56h, xx

xx: Highest possible value for ml/s + 1

4.5 Messages / Replies to Commands

The messages / replies are sent by the station in response to such frames which do not request a return value, or which have errors.

F1h Command understood
All is OK, command has been understood and run.

F2h Command unknown
Reply F2h, xx - command xx is undefined.

F3h Parameter is not permissible
The transmitted parameter is out of range.

F4h Syntax error
The frame is either too short or too long (missing characters, or too many).

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51h Station mode

Query: 51h
Reply: 51h, xx

xx: 0 = IDLE, standstill, idle
 1 = MANUAL, manual displacement
 2 = DISP, dispensing, station ejects
 3 = TOPD, station is in the top delay state
 4 = ASPI, Station is sucking in
 5 = PRED = Station is in the pre-delay state
 6 = USER, dialogue with the user is running

52h Number of repeats

Query: 52h
Reply: 52h, value

Value: Number of repeats not yet completed,
 8 digit BCD value (4 bytes)

53h End of last profile

Query: 53h
Reply: 53h, xx

xx: 0 = NORM, profile terminated normally
 1 = STOP, termination by operating the EMERGENCY STOP push-button
 2 = COMM, end by STOP command
 3 = PROG, end because of program error
 255 = RUN, profile is still running

44h LOCK Command

Command: 44h
Reply: F1h or F6h or F7h or F8h

Sets (if successful) the slave to the LOCKed state. In this state the station can not initiate a manual displacement or a user dialogue. The station can only be controlled via the interface. This command should be sent ahead of complex operations via the interface. An OPTIMAT® controller sends the command before processing a code list.

45h UNLOCK Command

Command: 45h
Reply: F1h

Cancels the LOCKed state of the slave. An OPTIMAT® controller sends the command after processing a code list.

46h Initiate BREAK

Command: 46h
Reply: F1h

Causes the corresponding slave to activate its EXBREAK\ line. This causes a BREAK in all connected slaves (including the one which was addressed). Therefore this command needs only to be sent to one slave in the network.

4Fh Test Command

Command: 47h
Reply: F1h

Causes all LEDs of the addressed slave to light up briefly; will not cause any operations or effects. This command may be used to check the link when establishing a link between computer and OPTIMAT®. Will possibly be no longer implemented in subsequent software releases.

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Example: Slave 15 is to be set to a dispensing volume of 1234.56 ml.

Frame:

8Fh, 56h, 34h, 12h, 00h, 0Dh

Frame end

Dispensing volume, 8 digits BCD,
least significant digit first

Slave address, binary OR 128

4.3 Sequencing Commands

Starting / stopping of processes

41h STOP

Command: 41h

Reply: F1h

Stops the current dispensing process

42h START

Command: 42h

Reply: F1h or F5h or F6h or F7h or F8h

Starts the dispensing process in the current memory state.

43h START (suction only)

Command: 43h

Reply: F1h or F5h or F6h or F7h or F8h

Starts the dispensing process in the current memory state; suction only, the dispensing process is cancelled after the top delay has elapsed.

21h	Set/query volume Value: Volume µl - from 0 to 9999999
22h	Set/query velocity Value: Velocity in 1/100 ml/s - from 0 to 9999
23h	Set/query velocity down Value: Velocity in 1/100 ml/s - from 0 to 9999
24h	Set/query pre-delay Value: Pre-delay in s - from 0 to 9999
25h	Set/query top delay Value: Top delay in s - from 0 to 9999
26h	Set/query repeat Value: Number of repeats - from 0 to 9999
27h	Set/query I/O type Value: Input = Output = 00 00 00 00h Input > Output = 00 00 00 01h
28h	Set/query acceleration up Value: Acceleration - from 0 to 999
29h	Set/query acceleration down Value: Acceleration - from 0 to 999
2Ah	Set/query deceleration up Value: Deceleration - from 0 to 999
2Bh	Set/query deceleration down Value: Deceleration - from 0 to 999
2Ch	Set/query next memory Value: Next memory - from 0 to 15

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To date, the following system parameters have been defined:

- 10: Sound signal during delay phases
off = 0, on = 1, default : 1
- 80: Signal lead for the solenoid valves
range: 0 ... 200, default: 0
- 81: Speed when driven manually
range: 0 ... 30, default: 20
- 90: CAN identifier
range: 0 ... 255, default 230
- 91: CAN bus timing 1
range: 0 ... 255, default: 65
- 91: CAN bus timing 2
range: 0 ... 255, default: 28

4.2 Memory Related Commands / Queries

Manipulation / analysis of the currently set memory.

All dialogue with 2x commands is run as follows:

Query: 2xh
Reply: 2xh, value

Set: 2xh, value
Reply: F1h

Value: Parameter value - from/to, see individual commands

All values are transmitted by 8 digit BCD numbers through 4 bytes (given in the diagram only as 'value'). The value is transmitted with the least significant digits first. Here the reply of an error-free dialogue is stated as the reply. For replies in the case of errors, see next chapter.

4 Commands of the Protocol

The numbers of the commands are indicated in the following in hexadecimal form. The corresponding ASCII character is indicated after the hex. value.

4.1 Station Related Commands/Queries

11h Select/query memory

Query: 11h

Reply: 11h, xx

Set: 11h, xx

Reply: F1h

xx: memory number - 01h ... 15h (1 ... 15 BCD)

12h Select/query size of the pump

Query: 12h

Reply: 12h, xx

Set: 12h, xx

Reply: F1h if all is OK
 F9h when the automatic system has detected a
 different flask

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
3.3 Addressing

If a station receives a frame which contains its own station address (see below), it will process that frame and will return a frame containing its own station number. You must under all circumstances avoid the assignment of the same number to two stations.

As a rule the following applies (as seen from the master):

For Transmitting:

Station address is the number of the slave for which the command/query is intended.

For receiving:

Station address is the number of the slave which has sent the response (must be the same number as used in the preceding command / query frame).

The station addresses are binary coded and are transmitted with bit 7 set (logic OR operation with number 128).

3.4 Structure of the Frames

Frame end character (Cr, 0xD)

Frame content

Station address

Each command is sent in a frame having 2 to n characters. At RSxxx the first character is the station address and a CR (<U, 0Dh) is sent as the frame end character.

All values are transmitted in the BCD format. Here always all digits must be stated (e.g. '00000000' instead of only '0').

In the case of CAN, the entire (RSxxx-) frame is sent in the data field, only the end character is omitted. A value which still remains to be defined is used as the identifier for all stations, i.e. each station receives all frames.

All pins marked “reserved” have been reserved for future extensions, but these pins are not yet being used. No external signals must be connected to these.

Only the signals of pins 2 and 3 are used for communication purposes.

2.2 Interface Mode

	U	C	S
none	yes	no	no
RS 232	yes	yes ¹⁾	yes ¹⁾
RS 485	yes	yes	yes
CAN	yes	yes	yes

¹⁾ with one slave only

3 Protocol

3.1 Processing / Transmission Rights

The protocol used is designed as a master-slave protocol, i.e. a station - the master (a PC or an OPTIMAT® in the controller mode) - sends a frame to a slave and waits for a response from the slave. Only a slave which previously has received a frame from the master may send one frame back to the master. In the master itself a routine must be implemented for the purpose of detecting a time-out.

3.2 Station Numbers

In order to discern the individual slaves in a network, a station number is introduced (range 1 ... 99). The number of a station is set up via the keyboard of each station. Each number may only be assigned to a slave one. Errors in assigning the station number can not be detected by the current software release and will cause unpredictable effects when processing commands via the interface. The numbers need not be assigned consecutively.

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All pins marked "reserved" have been reserved for future extensions, but these pins are not yet being used. No external signals must be connected to these.

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2.1.3 CAN (Controller Area Network)

As the CAN controller a PCA82C200 standalone controller (Basic CAN, CAN spec. 1.2) is used. The driver is a modified RS 485 driver. The bus timing and the identifiers used are defined through system parameters 90-92.

90	- Identifier	(default: 230)
91	- Bus timing 1	(default: 65)
92	- Bus timing 2	(default: 28)

A female 25-way sub-D connector is used.

Pin assignment:

2	CAN Hi	Line "A" of the RS 485 driver
3	CAN Lo	Line "B" of the RS 485 driver
7	GND	- signal ground
9	+5 V, 50 mA	
10	reserved	
11	MSTOP \	- stop motor (STOP push-button)
12	MLED_L	- Motor moves to the left
13	MLED_R	- Motor moves to the right
14	reserved	
15	reserved	
16	reserved	
17	reserved	
18	EXBREAK\	- external BREAK signal
23	reserved	
24	MHSL\	- Motor moves - manually operated - to the left
25	MHS_R\	- Motor moves - manually operated - to the right

Pin numbers which are not listed have not been assigned within the instrument and are not connected.

Only the signals of pins 2, 3 and 7 are used for communication purposes. The data are transmitted at a signal level of ± 10 V.

With a RS 232 interface you may only arrange point-to-point links, i.e. you may only control one OPTIMAT® via a PC interface or an other OPTIMAT®.

2.1.2 RS 485

The serial interface in accordance with RS 485 uses two lines A and B for communication purposes.

The transmission parameters corresponds to those of the RS 232 interface.

A female 25-way sub-D connector is used.

Pin assignment:

2	Line "A" of the RS 485 channel
3	Line "B" of the RS 485 channel
7	GND - signal ground
9	+5 V, 50 mA
10	reserved
11	MSTOP \ - stop motor (STOP push-button)
12	MLED_L - Motor moves to the left
13	MLED_R - Motor moves to the right
14	reserved
15	reserved
16	reserved
17	reserved
18	EXBREAK\ - external BREAK signal
23	reserved
24	MHSL\ - Motors moves - manually operated - to the left
25	MHS_R\ - Motors moves - manually operated - to the right

Pin numbers which are not listed have not been assigned within the instrument and are not connected.

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2.1.1 RS 232C

The serial interface in accordance with RS 232C uses for communication purposes only the signals RxD, TxD and GND. Due to small amount of data, hardware handshaking or a low level software handshake (XON/XOFF) is not required.

Specifications of the interface:

- Transmission rate: 9600 Baud
- Data bits: 8
- Parity: none
- Stop bits: 1

A female 25-way sub-D connector is used which is normally wired as far as the pins are used at all.

Pin assignment:

- 2 TxD - transmitted data
- 3 RxD - received data
- 7 GND - signal ground
- 9 + 5 V, 50 mA
- 10 reserved
- 11 MSTOP \ - stop motor (STOP push-button)
- 12 MLED_L - Motor moves to the left
- 13 MLED_R - Motor moves to the right
- 14 reserved
- 15 reserved
- 16 reserved
- 17 reserved
- 18 EXBREAK\ - external BREAK signal
- 23 reserved
- 24 MHSL\ - Motor moves - manually operated - to the left
- 25 MHS_R\ - Motor moves - manually operated - to the right

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
1.1 Configurations without a PC

One station alone - UNIQUE


Two or more new stations (bus) also minimum version with 2 stations and RS 232

1.2 Configurations with a PC

PC + one or several new station(s) (bus). New command set, also minimum version with one station and RS 232

2 Interfaces

2.1 Types of Interface

The OPTIMAT may be equipped with 3 different interfaces:

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yy: Parameterwert, 0 ... 255, BCD - codiert, 2 Bytes

ANNEX A

Maximum speed for intake and dosing procedure

Pump size	V _{max} in ml/s
1 ml	0,55
2 ml	1,90
10 ml	3,35
30 ml	9,60
50 ml	15,20
100 ml	24,30
200 ml	31,30
300 ml	31,80

Attention:

When entering invalid values faster than maximum speed, no error message will be displayed! The dosing station will automatically set the value of the allowed maximum speed. Beware of setting the speed index to zero: the station would not work.

ANNEX B

FORTUNA® OPTIMAT® 2 - Interface Protocol

1 Configurations / Types of Station

U = UNIQUE

No or any kind of interface
Interface is not used

C =

Controller on the bus
RS 485 or CAN interface (RS 232 with restrictions)
Sends a new command set to new stations

S =

Slave on the bus
RS 485 or CAN interface (RS 232 with restrictions)
Only together with controller or PC (new software)

8 a) Maintenance

By its mainly closed construction the dosing station FORTUNA® OPTIMAT® 2 is protected from contaminations. All mechanical parts used are nearly maintenance-free. In order to guarantee a long-lasting functioning, we recommend a regular technical control either by yourself or by the manufacturer in a 2 years period.

8 b) Sterilization

Removal of dosing pumps (see manual page 34):

Remove the screw of the upper pump holder and remove the pump from the lower pump holder by opening the PTFE screw of the pump.

Dismantling of dosing pumps: Pull the dosing piston out of the dosing cylinder.

Sterilization: The dismantled dosing pump can be sterilized at +121° C (+248°F)/1 bar (max.).

Cooling: After the sterilization the dismantled dosing pump should be allowed to cool down very slowly at room temperatur (minimum 1 hour).

Assembly of dosing pumps: After cooling the dosing piston may be inserted into the dosing cylinder.

Installation of dosing pumps: Please follow the instructions on page 6 and 7.

9 Warranty

The manufacturer guarantees each dosing station for a period of 2 years excluding all claims for damage, conversion or other claims. The dosing station FORTUNA® OPTIMAT® 2 is maintenance-free and extremely reliable. However, if any trouble arises, the manufacturer reserves the right to repair or exchange faulty parts. The manufacturer cannot be held responsible for damages or troubles caused by improper handling or maintenance. All claims should be reported to the supplier. The warranty period shall not be extended as a result of repair under warranty, neither for the replacement parts nor for the dosing station itself. The warranty shall become null and void, if repairs have been carried out by unauthorized persons.

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5 Troubleshooting - Possible trouble and resulting actions

Trouble:	No message on the display after switching on the dosing station.
Action:	Check whether station shows damage due to transportation and whether connection is correct and check the fuse. If there's no damage visible and/or the cable is connected, please contact your dealer or technical service.
Trouble:	The dosing pump dispenses too much or not enough.
Action:	Check pump installation. Check dosing program. Check tightness of connection tubes. Check correct function of dosing pump. Run the motor manually down with the motor keys.
Trouble:	The dosing station does not start.
Action:	Check program. Avoid zero inputs. Switch off station and switch on again after 5 minutes. If station does not react, contact technical service.
Trouble:	The dosing station counts or works wrongly.
Action:	Interrupt the program by pressing the emergency STOP key. Press F2 SETUP and select F3 PUMP. Follow the instructions on the display, check size of the pump and confirm all input with F4 OK or F4 ENTER. Check the pump. By pressing F3 EXIT the installation procedure can be quitted. After a new pump setup check programming again.

6 Titration

A titration can be made by defining the operation speed combined with the delay times.

7 Storage of program

Program memories ensure that all values remain stored also when station is switched off or in case of power failure. However, it is recommended to check the programs when the station is not in use or switched off the line for several weeks or after mains failure.

4.4 Memory selection

To modify or create a program for unit 0 (controller) follow the steps as described on page 8 ff. To modify or create a program for unit >00 (slaves) first do step 4.3. **Memory selection**, to choose the number of the unit and follow the steps as described on page 8 ff.

4.5 Starting the selected dosing program

Press F1 EXECUTE to start the selected dosing program.
Press the Motor STOP key to interrupt the dosing procedure or in an emergency.
After the dosing procedure has been finished, the display shows the start mode, functions F1 - F4 may be selected.

4.6 System edit mode for double station programming (only for technical service or adviced operators!)

If the FORTUNA® OPTIMAT® 2 was ordered as a double station this programming is not necessary and already done at the manufacturer. If the programming was changed by pressing F4 EDIT the double station function might become invalid!

Before entering the system edit mode please start with
11.1 Defining the status of the dosing station.

Press F4 EDIT to enter the system programming mode. By pressing F4 EDIT again the cursor will go to the first program line (if there is any). By pressing F4 ENTER input can be confirmed or values can be changed by pressing the flashing edit mode keys.
Press F3 QUIT to close the menu.

If the double station shall be used as 2 single stations please do step 11.1 and change value for the dosing station to «UNIQUE».

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4 Operation as a DOUBLE station

To use the FORTUNA® OPTIMAT® 2 as a double station 2 units must be equipped with an interface and connected with a cable (option).

4.1 Defining the status of the dosing station

The first FORTUNA® OPTIMAT® 2 must be defined as a CONTROLLER, the second unit as a SLAVE. All further units in the chain will also be defined as SLAVE.

The display shows the 4 functions **F1 EXECUTE**, **F2 SETUP**, **F3 MEM+PRG** and **F4 EDIT**.

Press **F2 SETUP** to enter the controller/slave status.

Press **F2 STATUS**.

status: Contr. (Slave / Unique)	Change the value by pressing the flashing edit key and confirm with F4 ENTER .
---------------------------------	---

4.2 Setup of the pump

For the setup of the pump please see page 6 ff.

4.3 Memory selection

Press **F3 MEM+PRG** and follow the program, change values by pressing the flashing edit keys and confirm with **F4 ENTER**.

sta: XX (Station) Input of the number of the dosing station.
(The controller must always be no. 0, all further slaves are numbered from 01 - 99)

mem: XX Select your program number (1-15).

To quit the programming press **F3 EXIT**.

Warning: By pressing **F3 ENTER** the programming might be not complete!

acceler. ↓: XXX	(Acceleration) Input of a value between 0 (no ramp) and 800 (slow ramp) to create a speed profile (ramp) when the dispensing procedure starts from the top downwards.
deceler. ↑: XXX	(Deceleration) Input of a value between 0 (no ramp) and 800 (slow ramp) to create a speed profile (ramp) when the dispensing procedure comes to the top.
deceler. ↓: XXX	(Deceleration) Input of a value between 0 (no ramp) and 800 (slow ramp) to create a speed profile (ramp) when the dispensing procedure ends downwards.
next memory: XX	Input of a next memory to create a dispensing chain (connecting jobs). If the input is zero no program will be connected and the job ends.

(End of programming.)

3.5 Starting the selected dosing program

Press **F1 START** to start the selected dosing program.
Press the Motor **STOP** key to interrupt the dosing procedure or in an emergency.
After the dosing procedure has been finished, the display shows the start mode, functions **F1 - F4** may be selected.

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predelay: XXXXs

1 repetition: Input of the delay time between pressing **F1 START** and the first dispensing.

>1 repetition: Input of the delay time between pressing **F1 START** and the first dispensing and between each dispensing repetition.

topdelay: XXXXs

Input of the time to delay the dispensing procedure when the piston of the pump is in the top position.

(This feature is useful to ensure that the valve ball closes well when liquids with high viscosity are dispensed.)

repeat: XXXXXx

Input of repetitions.

Attention: No function when zero!

typ: Input=Output

Input of volume (vol) is equal to the dispensed volume.

Input>Output

The pump sucks always the nominal volume when it is empty or there is not enough liquid left in the pump to do the next dispensing.

The volume dispensed is referred to the volume input (vol).

Modify this function by pressing the  key.

acceler. : XXX

(Acceleration)

Input of a value between 0 (no ramp) and 800 (slow ramp) to create a speed profile (ramp) when the dispensing procedure starts upwards.

3.3 Memory selection

Press **F3 MEMORY** and select your program number (1-15) by using the edit keys **C**. Confirm your selection with **F4 ENTER**. Press **F3 EXIT** to abort.

3.4 Setting up the program of a memory

Press **F4 PROGRAM** and follow the program steps. Change the given values by means of the edit keys **C**. Confirm all input with **F4 ENTER**. Press **F3 EXIT** to abort programming.
The flashing lights indicate the possible keys which can be used for input. Not flashing lights can not be used for input.

ATTENTION: By pressing **F3 EXIT** all previous input confirmed by **F4 ENTER** remains stored except the last data not confirmed by **F4 ENTER**.

Explanation of the programming steps (function display):

vol:	XXXX.XXXml	(Volume) Input of dosing volume.
velo. ↑:	XXX.XXXml/s	(Velocity upwards in ml per second) Input of the intake speed.
velo. ↓:	XXX.XXXml/s	(Velocity downwards in ml per second) Input of the discharge speed.

Table of minimum/maximum speed
please find in annex A.

Attention: If the speed value input is
out of range, the speed
rate will be set
automatically to the max.

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Fix piston

Now fix the head of the piston by tightening the screw. Press **F4 OK**.

The standard display with four functions is shown.

Connect the intake and discharge tube to the valve block of the pump.

Option:

If an electronic handprobe is used connect also the cable of the handprobe to the resp. interface on the back of the dosing station.

Now the dosing station is ready to start or the memories can be programmed.

3 Operation as a UNIQUE (single) station

3.1 Connection

Please check whether delivery is complete according to the delivery note. Connect the dosing station FORTUNA® OPTIMAT® 2 to the power supply system by means of the connection cable. Use only original FORTUNA® parts. The station recognizes interfaces automatically and whether the voltage is 220 or 110 V and sends the message 'OPTIMAT® PRG-VERSION X.XX' as well as several short sound signals.

3.2 Installation of pump in SETUP

Press **F2 SETUP** and select **F3 PUMP**. Follow the instructions on the display as shown below. Confirm with **F4 OK**. Press **F3 EXIT** to quitt the menu.

GB

- | | |
|---------------------|--|
| Install new pump | Press F4 OK to confirm. |
| Remove old pump | Remove the screw at F and remove the pump from G by opening the PTFE screw of the pump. Press F4 OK . |
| Fix base of pump | Fix base of new pump at G. Do not yet fix head of pump. Press F4 OK . |
| Pump : XXX ml | Press ↑ or ↓ to set the size of the new pump. Press F4 ENTER . |
| Move piston down | by hand to ensure that the piston is at the zero position. Press F3 READY . |
| Calibrating standby | calibrates the motor head automatically to the zero position. |

Preface

We thank you for the purchase of the FORTUNA® OPTIMAT® dosing station. You did a good choice. With more than 20 years of experience in making dosing systems we guarantee a high standard of product quality and safety.

To guarantee quality and safety of the dosing station also for a long term, a correct operating and maintenance is necessary.

Therefore please read this operating manual carefully before using the dosing station and take notice of the warning instructions as well as the corresponding regulations.

The delivery contains: 1 FORTUNA® OPTIMAT® 2 dosing station with integrated programming terminal.
Please see our catalog for the options.

1 The world of FORTUNA® OPTIMAT® 2

150.000-2	FORTUNA® OPTIMAT® 2 COMFORT with integrated programming terminal
150.000-3	FORTUNA® OPTIMAT® 2 COMFORT / S with integrated programming terminal and connection port RS 232 (standard), RS 485 or CAN
Optional:	
150.000-10	Interface RS 232 (only for cat. no. 150.000-2)
150.000-20	Interface RS 485 (only for cat. no. 150.000-2)
150.000-30	Interface CAN (only for cat. no. 150.000-2)

The dosing station FORTUNA® OPTIMAT® 2 is suitable for one OPTIMAT® pump only (please see our catalog for accessories). The dosing station 150.000-3 OPTIMAT® 2 COMFORT / S can be networked with further dosing stations OPTIMAT® 2 COMFORT / S or other hardware (e. g. printer, balance) via a cable (option).

2 Programming Terminal

A Function keys

The function keys F1 - F4 change according to the program.

B LCD-display

The display has two lines with four fields for the function key commands and is illuminated.

C Cursor keys

By pressing ↑ and ↓ the given value can be increased resp. decreased.
By pressing ← and → the cursor can be moved to the resp. digits.

D Memory LED

15 LED to indicate the selected program.

E Motor control

Interrupts current program.

The motor can be manually moved by the ↑ and ↓ keys.

F Pump holder (head)

The head holder of the pump is connected to the motor unit and the piston of the pump.

G Pump holder (base)

The base holder fixes the valve block of the pump.